

Serum Drug Screen Comprehensive (SDSC) Drug List

DRUGS	HPLC			
	Immunoassay	Basic	Polar (#)	GC/FID (#)
Acetaminophen			#	
Alprazolam				#
Amiodarone		x, ϕ		
Amitriptyline		#		
Amphetamine	x			
Barbiturates**	x			
Benzodiazepines*	x			#
Bupropion (metabolite)		x		
Caffeine			#	
Carbamazepine			#	
Chlordiazepoxide				#
Chlorpheniramine		x		
Chlorpromazine		x		
Citalopram		x		
Clomipramine		x		
Clonazepam				#
Clozapine		x		ϕ
Cocaine Metabolite	x			
Codeine				ϕ
Cyclobenzaprine		x		
Desipramine		#		
Diazepam				#
Diltiazem		x		
Diphenhydramine		x		
Doxepin		#		
Duloxetine		x		
Ethanol	x			#
Felbamate			#	
Fluoxetine		x, ϕ		
Hydrocodone				ϕ
Hydromorphone				ϕ
Imipramine		#		
Lamotrigine				ϕ
Levetiracetam				ϕ
Lorazepam				#
Maprotiline		x		
MDA	x			
MDMA	x			
Meperidine		x		
Mephobarb			#	
Mesoridazine		x		
Methadone	x			ϕ
Methamphetamine	x			
Metoprolol		x		

DRUGS	HPLC			
	Immunoassay	Basic	Polar (#)	GC/FID (#)
Morphine				ϕ
Mirtazapine		x		
N-Desalkylflurazepam				#
Nordiazepam				#
Nortriptyline		#		
Olanzapine		x		
Opiates	x, §			ϕ
Oxazepam				#
Oxcarbazepine			#	
Oxycodone				ϕ
Oxymorphone				ϕ
Pentazocine		x		
Pentobarbital			#	
Phencyclidine	x			
Phenobarbital			#	
Phenytoin			#	
Promazine		x		
Promethazine		x		
Propafenone		x		
Propoxyphene	x	x		
Propranolol		x		
Protriptyline		x		
Quinine/Quinidine		x		
Salicylates			#	
Sertraline		x		
Temazepam				#
THC Metabolites	x			
Thioridazine		x		
Tramadol		x, ϕ		
Trazodone				ϕ
Trifluoperazine		x		
Trimipramine		x		
Venlafaxine		x		
Verapamil		x		

KEY
x = Qualitative Results Reported
= Quantitative Results Reported
ϕ = Quantitation available upon request at additional charge. Contact Lab.
§ = Opiates detected only at levels over 100 ng/mL. Both free & conjugated.
*Most benzodiazepines detected. Most likely no Flunitrazepam or triazolam.
**Quantitated Barbiturates are phenobarbital and pentobarbital.

