

Update Notes

Update Summary

New Test Activation	4/23/2019	BNUGG - "Brazil Nut Food IgG"
New Test Activation	4/23/2019	CNNGG - "Cinnamon IgG"
New Test Activation	4/23/2019	COCGG - "Coconut IgG"
New Test Activation	4/23/2019	JACHI - "Jalapeno/Chipotle IgE"
New Test Activation	4/23/2019	STPII - "Staphylococcal Panel II"
New Test Activation	4/23/2019	WAYEP - "Wall-Eyed Pike IgE"
New Website Listing	4/23/2019	ILFG1 - "Insulin-Like GF1 Binding Protein"
New Website Listing	4/22/2019	ULIP - "Lipase-Urine, Random"
Update Existing Test	4/8/2019	CHRMS - "Chromium, Serum"
Update Existing Test	4/15/2019	PAI1G - "Plas Act Inhib (PAI-1) Genotype"
Update Existing Test	4/15/2019	PHGAM - "Phosphatidylserine Antibodies (IgG IgM IgA)"
Update Existing Test	4/8/2019	UAMGF - "Amphetamine, Urine, Forensic GC/MS Confirmation (Quantitative)"
Update Existing Test	4/8/2019	UBAGF - "Barbiturate, Urine, Forensic, GC/MS Confirmation (Quantitative)"
Update Existing Test	4/8/2019	UBZGF - "Benzodiazepine, Urine, Forensic Confirmation (Quantitative)"
Update Existing Test	4/8/2019	UCAMP - "Clin Urine Amphetamine Confirm"
Update Existing Test	4/8/2019	UCBAR - "Clin Urine Barbiturate Confirm"
Update Existing Test	4/8/2019	UCBEN - "Clin Urine Benzodiazepine Confirm"
Update Existing Test	4/8/2019	UCBUP - "Clin Urine Buprenorphine Confirm"
Update Existing Test	4/8/2019	UCCOC - "Clin Urine Cocaine Confirm"
Update Existing Test	4/8/2019	UCETH - "Clin Urine Ethanol Confirm"
Update Existing Test	4/8/2019	UCFNT - "Clin Urine Fentanyl Confirm"
Update Existing Test	4/8/2019	UCMPE - "Clin Urine Meperidine Confirm"
Update Existing Test	4/8/2019	UCMTH - "Clin Urine Methadone Confirm"
Update Existing Test	4/8/2019	UCOGF - "Cocaine, Urine, Forensic GC/MS Confirmation (Quantitative)"
Update Existing Test	4/8/2019	UCOPT - "Clin Urine Opiate Confirm"
Update Existing Test	4/8/2019	UCPCP - "Clin Urine Phencyclidine Confirm"
Update Existing Test	4/8/2019	UCPRP - "Clin Urine Propoxyphene Confirm"

Update Existing Test	4/8/2019	UCTRM - "Clin Urine Tramadol Confirm"
Update Existing Test	4/8/2019	UDS01 - "Drug Screen, Urine Comprehensive"
Update Existing Test	4/8/2019	UETG3 - "EtG Screen w/ EtG/EtS Confirmation"
Update Existing Test	4/8/2019	UETHF - "Ethanol, Urine, Forensic Confirmation (Qualitative)"
Update Existing Test	4/8/2019	UFEN5 - "Clin Urine Fentanyl Scrn w/reflex Conf"
Update Existing Test	4/8/2019	UMAC1 - "6-Acetylmorphine, Clinical, Quantitative"
Update Existing Test	4/8/2019	UMAF1 - "6-Acetylmorphine, Forensic, Quantitative"
Update Existing Test	4/8/2019	UMEP5 - "Clin Urine Meperidine Screen w/reflex Conf"
Update Existing Test	4/8/2019	UMTGF - "Methadone, Urine, Forensic, GC/MS Confirmation"
Update Existing Test	4/8/2019	UOPGF - "Opiate, Urine, Forensic, LC/MS/MS Confirmation (Quantitative)"
Update Existing Test	4/8/2019	UPCGF - "Phencyclidine, Urine, Forensic, GC/MS Confirmation (Quantitative)"
Update Existing Test	4/8/2019	UPRGF - "Propoxyphene, Urine Forensic GC/MS Confirmation (Quantitative)"
Update Existing Test	4/8/2019	UTRM5 - "Clin Urine Tramadol Scrn w/reflex Conf"

New Test Activation

Effective Date	4/23/2019
Name	Brazil Nut Food IgG
Code	BNUGG
CPT Code(s)	86001 RUO
Notes	

Specimen Requirements

Specimen Required	Please see the Allergen, IgG website listing for specimen requirements and other information.
-------------------	---------------------------------------------------------------------------------------------------------------

Performing Information

Performing Laboratory	Viracor Eurofins
-----------------------	------------------

Interface Information

Legacy Code ¹	BNUGG		
Interface Order Code	3300008		
Result Code	Name	LOINC Code	AOE/Prompt ²
3300008	Brazil Nut Food IgG	63081-4	No

New Test Activation			
Effective Date	4/23/2019		
Name	Cinnamon IgG		
Code	CNNGG		
CPT Code(s)	86001 RUO		
Notes			
Specimen Requirements			
Specimen Required	Please see the Allergen, IgG website listing for specimen requirements and other information.		
Performing Information			
Performing Laboratory	Viracor Eurofins		
Interface Information			
Legacy Code ¹	CNNGG		
Interface Order Code	3300009		
Result Code	Name	LOINC Code	AOE/Prompt ²
3300009	Cinnamon IgG	60367-0	No

New Test Activation			
Effective Date	4/23/2019		
Name	Coconut IgG		
Code	COCGG		
CPT Code(s)	86001 RUO		
Notes			
Specimen Requirements			
Specimen Required	Please see the Allergen, IgG website listing for specimen requirements and other information.		
Performing Information			
Performing Laboratory	Viracor Eurofins		
Interface Information			
Legacy Code ¹	COCGG		
Interface Order Code	3300018		
Result Code	Name	LOINC Code	AOE/Prompt ²
3300018	Coconut IgG	63084-8	No

New Test Activation			
Effective Date	4/23/2019		
Name	Jalapeno/Chipotle IgE		
Code	JACHI		
CPT Code(s)	86003		
Notes			
Specimen Requirements			
Specimen Required	Please see the Allergen, IgE website listing for specimen requirements and other information.		
Performing Information			
Performing Laboratory	Viracor Eurofins		
Interface Information			
Legacy Code¹	JACHI		
Interface Order Code	3300044		
Result Code	Name	LOINC Code	AOE/Prompt²
3300045	Jalapeno/Chipotle IgE	11189-8	No
3300046	Class	21453-6	No

New Test Activation			
Effective Date	4/23/2019		
Name	Staphylococcal Panel II		
Code	STPII		
CPT Code(s)	86003 x 2, 82785		
Notes			
Specimen Requirements			
Specimen Required	Draw blood in a SST. Centrifuge, separate serum from cells, and send 1.5 mL serum (1.0 mL minimum) room temperature in a screw-capped plastic vial.		
Alternate Specimen	Serum: Red-top		
Rejection Criteria	Gross lipemia		
Stability	Room temperature: 8 days; Refrigerated: 8 days; Frozen: 4 weeks		
Performing Information			
Methodology	ImmunoCAP® FEIA		
Reference Range	By report		
Performed Days	Varies		
Turnaround Time	4-7 days		
Performing Laboratory	Viracor Eurofins		
Interface Information			
Legacy Code ¹	STPII		
Interface Order Code	3300051		
Result Code	Name	LOINC Code	AOE/Prompt ²
3300052	Total IgE	19113-0	No
3300053	Staphylococcal Enterotoxin A IgE	25821-0	No
3300054	Staphylococcal Enterotoxin B IgE	25822-8	No

New Test Activation			
Effective Date	4/23/2019		
Name	Wall-Eyed Pike IgE		
Code	WAYEP		
CPT Code(s)	86003		
Notes			
Specimen Requirements			
Specimen Required	Please see the Allergen, IgE website listing for specimen requirements and other information.		
Performing Information			
Performing Laboratory	Viracor Eurofins		
Interface Information			
Legacy Code ¹	WAYEP		
Interface Order Code	3300031		
Result Code	Name	LOINC Code	AOE/Prompt ²
3300032	Wall-Eyed Pike IgE	7761-0	No
3300033	Class	21460-1	No

New Website Listing			
Effective Date	4/23/2019		
Name	Insulin-Like GF1 Binding Protein		
Code	ILFG1		
CPT Code(s)	83519		
Notes			
Specimen Requirements			
Specimen Required	Overnight fasting is preferred. Draw blood in a SST. Centrifuge, remove serum from cells and send 1.0 mL serum (0.2 mL minimum) refrigerated in a screw-capped plastic vial.		
Rejection Criteria	Gross hemolysis, gross lipemia		
Stability	Room temperature: 24 hours; Refrigerated: 7 days; Frozen: 28 days		
Performing Information			
Methodology	Radioimmunoassay		
Reference Range	Adults: 5-34 ng/mL Pediatric: 5-9 years: 15-95 ng/mL 10-14 years: 8-64 ng/mL 15-18 years: 5-40 ng/mL		
Performed Days	Tuesday		
Turnaround Time	4-10 days		
Performing Laboratory	Quest SJC		
Interface Information			
Legacy Code¹	ILFG1BP		
Interface Order Code	3503750		
Result Code	Name	LOINC Code	AOE/Prompt²
3503750	Insulin-Like GF1 Binding Protein	12722-5	No

New Website Listing			
Effective Date	4/22/2019		
Name	Lipase-Urine, Random		
Code	ULIP		
CPT Code(s)	83690		
Notes			
Specimen Requirements			
Specimen Required	Send 2.0 mL random urine (1.0 mL minimum) refrigerated in a screw-capped plastic vial.		
Alternate Specimen	Urine: aliquot from 24 hour collection.		
Stability	Room temperature: 5 days; Refrigerated: 5 days; Frozen: 21 days		
Performing Information			
Methodology	Spectrophotometry (SP)		
Reference Range	<=4 U/L		
Performed Days	Sunday, Tuesday-Friday		
Turnaround Time	4-6 days		
Performing Laboratory	Quest SJC		
Interface Information			
Legacy Code ¹	ULIP		
Interface Order Code	3504380		
Result Code	Name	LOINC Code	AOE/Prompt ²
3504380	Lipase-Urine	10888-6	No

Update Existing Test	
Effective Date	4/8/2019
Name	Chromium, Serum
Code	CHRMS
Interface Order Code	3422860
Legacy Code	CHRMQ
Notes	
Required Testing Changes	
Specimen Required	<p>Patient preparation: Patient should refrain from taking vitamins, mineral and herbal supplements at least 1 week prior to specimen collection.</p> <p>Draw one dark blue no additive tube and discard. Draw second dark blue no additive tube, allow to clot in an upright position. Spin and send 2.0 mL serum (1.0 mL minimum) refrigerated in white capped trace metal free plastic screw-capped vial.</p> <p>Tests performed on a specimen submitted in a nontrace element tube or nonacid-washed/nonmetal-free container may not accurately reflect the patients level. If a nontrace element tube/container is received, it will be accepted for testing; however, elevated results are reported with the message that resubmission with a trace element tube/container is recommended.</p>

Update Existing Test	
Effective Date	4/15/2019
Name	Plas Act Inhib (PAI-1) Genotype
Code	PAI1G
Interface Order Code	3426020
Legacy Code	PAI1GQ
Notes	
Required Testing Changes	
Alternate Specimen	Whole blood: Sodium or lithium heparin, ACD A, ACD B, dark blue EDTA
Methodology	Polymerase Chain Reaction and Detection

Update Existing Test

Effective Date	4/15/2019
Name	Phosphatidylserine Antibodies (IgG IgM IgA)
Code	PHGAM
Interface Order Code	3703180
Legacy Code	PHOSGMASP
Notes	

Required Testing Changes

Specimen Required	Draw blood in a light blue 3.2% sodium citrate tube. See appendices for coagulation test collection instructions. Send 1.0 mL platelet - poor plasma (0.5 mL minimum) refrigerated in a screw-capped plastic vial.
Alternate Specimen	Serum: Red-top or SST

Update Existing Test

Effective Date	4/8/2019
Name	Amphetamine, Urine, Forensic GC/MS Confirmation (Quantitative)
Code	UAMGF
Notes	

Required Testing Changes

Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.
--------------------------	--------------------------------------------------------------------------------------------------------------------------------

Update Existing Test

Effective Date	4/8/2019
Name	Barbiturate, Urine, Forensic, GC/MS Confirmation (Quantitative)
Code	UBAGF
Notes	

Required Testing Changes

Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.
--------------------------	--------------------------------------------------------------------------------------------------------------------------------

Update Existing Test	
Effective Date	4/8/2019
Name	Benzodiazepine, Urine, Forensic Confirmation (Quantitative)
Code	UBZGF
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send in 2.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Amphetamine Confirm
Code	UCAMP
Interface Order Code	1846150
Legacy Code	UCAMP
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Barbiturate Confirm
Code	UCBAR
Interface Order Code	1846500
Legacy Code	UCBAR
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Benzodiazepine Confirm
Code	UCBEN
Interface Order Code	1846800
Legacy Code	UCBEN
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Buprenorphine Confirm
Code	UCBUP
Interface Order Code	1846300
Legacy Code	UCBUP
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Cocaine Confirm
Code	UCCOC
Interface Order Code	1847100
Legacy Code	UCCOC
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Ethanol Confirm
Code	UCETH
Interface Order Code	1847200
Legacy Code	UCETH
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (2.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Fentanyl Confirm
Code	UCFNT
Interface Order Code	3000017
Legacy Code	UCFNT
Notes	
Required Testing Changes	
Specimen Required	Send 10.0 mL random urine (3.0 mL minimum) refrigerated in a screw-capped plastic container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Meperidine Confirm
Code	UCMPE
Interface Order Code	3000027
Legacy Code	UCMPE
Notes	
Required Testing Changes	
Specimen Required	Send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Methadone Confirm
Code	UCMTH
Interface Order Code	1847300
Legacy Code	UCMTH
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic container.

Update Existing Test	
Effective Date	4/8/2019
Name	Cocaine, Urine, Forensic GC/MS Confirmation (Quantitative)
Code	UCOGF
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Opiate Confirm
Code	UCOPT
Interface Order Code	1847600
Legacy Code	UCOPT
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.
Methodology	Liquid Chromatography - Tandem Mass Spectrometry

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Phencyclidine Confirm
Code	UCPCP
Interface Order Code	1847920
Legacy Code	UCPCP
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (8.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Propoxyphene Confirm
Code	UCPRP
Interface Order Code	1847840
Legacy Code	UCPRP
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Tramadol Confirm
Code	UCTRM
Interface Order Code	3000038
Legacy Code	UCTRM
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic container.

Update Existing Test	
Effective Date	4/8/2019
Name	Drug Screen, Urine Comprehensive
Code	UDS01
Interface Order Code	1820000
Legacy Code	UDS
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (15.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	EtG Screen w/ EtG/EtS Confirmation
Code	UETG3
Interface Order Code	1825390
Legacy Code	UETG3
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (2.0 mL minimum) refrigerated in a screw-capped plastic urine cup.

Update Existing Test	
Effective Date	4/8/2019
Name	Ethanol, Urine, Forensic Confirmation (Qualitative)
Code	UETHF
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (2.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Fentanyl Scrn w/reflex Conf
Code	UFEN5
Interface Order Code	3000011
Legacy Code	UFEN5
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic vial.

Update Existing Test	
Effective Date	4/8/2019
Name	6-Acetylmorphine, Clinical, Quantitative
Code	UMAC1
Interface Order Code	1836900
Legacy Code	UMAMCL
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Send 20.0 mL urine (8.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	6-Acetylmorphine, Forensic, Quantitative
Code	UMAF1
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Send 20.0 mL urine (8.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Meperidine Screen w/reflex Conf
Code	UMEP5
Interface Order Code	3000022
Legacy Code	UMEP5
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Methadone, Urine, Forensic, GC/MS Confirmation
Code	UMTGF
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Opiate, Urine, Forensic, LC/MS/MS Confirmation (Quantitative)
Code	UOPGF
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Phencyclidine, Urine, Forensic, GC/MS Confirmation (Quantitative)
Code	UPCGF
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (8.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Propoxyphene, Urine Forensic GC/MS Confirmation (Quantitative)
Code	UPRGF
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 20.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.

Update Existing Test	
Effective Date	4/8/2019
Name	Clin Urine Tramadol Scrn w/reflex Conf
Code	UTRM5
Interface Order Code	3000032
Legacy Code	UTRM5
Notes	
Required Testing Changes	
Specimen Required	Collect random urine. Mix well and send 10.0 mL urine (3.0 mL minimum) refrigerated in a screw-capped plastic urine container.