

Update Notes	

Update Summary		
Update Existing Test	3/28/2023	<a href="#">ACAB - "Centromere Antibodies"</a>
Update Existing Test	3/28/2023	<a href="#">NCA - "Neutrophil Cytoplasmic Ab"</a>
Update Existing Test	3/28/2023	<a href="#">RLTE4 - "Leukotriene E4, Random, Urine"</a>
Update Existing Test	3/28/2023	<a href="#">TLTE4 - "Leukotriene E4, 24 Hour, Urine"</a>
Inactivate Test Without Replacement	3/28/2023	<a href="#">CHGMS - "Chlamydia Serology, Serum"</a>
Inactivate Test Without Replacement	5/1/2023	<a href="#">PAIGM - "Parainfluenza Abs"</a>

Update Existing Test	
Effective Date	3/28/2023
Name	Centromere Antibodies
Code	ACAB
Interface Order Code	3006200
Legacy Code	ACAB
Notes	Update to alternate specimen, stability, methodology, reference range and performed days.
Required Testing Changes	
Alternate Specimen	Plasma: Lavender EDTA
Rejection Criteria	Heparinized plasma; hemolysis; lipemia
Stability	Room temperature: 8 hours Refrigerated: 14 days Frozen: Undetermined
Methodology	Enzyme Linked Fluorescent Immunoassay
Reference Range	<7 U/mL Negative 7-10 U/mL Equivocal >10 U/mL Positive
Performed Days	Monday - Friday

Update Existing Test	
Effective Date	3/28/2023
Name	Neutrophil Cytoplasmic Ab
Code	NCA
Interface Order Code	3007701
Legacy Code	NCA
Notes	Update to stability.
Required Testing Changes	
Stability	Room temperature: 8 hours <b>Refrigerated: 14 days</b> Frozen: 14 days

Update Existing Test	
Effective Date	3/28/2023
Name	Leukotriene E4, Random, Urine
Code	RLTE4
Interface Order Code	3800268
Legacy Code	
Notes	Update to stability.
Required Testing Changes	
Stability	Room temperature: 24 hours Refrigerated: 7 days <b>Frozen: 28 days</b>

Update Existing Test	
Effective Date	3/28/2023
Name	Leukotriene E4, 24 Hour, Urine
Code	TLTE4
Interface Order Code	3800271
Legacy Code	
Notes	Update to stability.
Required Testing Changes	
Stability	Room temperature: 24 hours Refrigerated: 7 days <b>Frozen: 28 days</b>

Inactivate Test Without Replacement	
<b>Effective Date</b>	3/28/2023
<b>Name</b>	Chlamydia Serology, Serum
<b>Code</b>	CHGMS
<b>Legacy Code</b>	CHGMAM
<b>Interface Code</b>	3800120
<b>Notes</b>	Suggested replacement(s): <a href="#">CHLAP – Chlamydia and Chlamydomphila Antibody Panel 3 IgG, IgA, IgM</a> <a href="#">CHRNA – Chlamydia Trachomatis Nucleic Acid by TMA</a>

Inactivate Test Without Replacement	
<b>Effective Date</b>	5/1/2023
<b>Name</b>	Parainfluenza Abs
<b>Code</b>	PAIGM
<b>Legacy Code</b>	PARA
<b>Interface Code</b>	3508400
<b>Notes</b>	Suggested replacement: <a href="#">PIPCR – Parainfluenza PCR</a>