

LABORATORY REPORT

Example Client, XYZ123 1234 Warde Road Ann Arbor MI 48108

EXAMPLE, REPORT W WX000003827 M 07/08/1978

Collected: 09/26/2023 14:24 Received: 09/26/2023 14:24

Drug Screen, Pain Manageme

<u>Test Name</u>	<u>Result</u>	<u>Flag</u>	<u>Ref-Ranges</u>	<u>Units</u>	<u>Site</u>
Amphetamine Screen	Negative				WMRL
Barbiturate Screen	Negative				WMRL
Benzodiazepine Screen	Negative				WMRL
Buprenorphine Screen	Negative				WMRL
Cocaine Screen	Negative				WMRL
Ethanol Screen	Negative				WMRL
Fentanyl	Negative				WMRL
Meperidine	Negative				WMRL
Methadone Screen	Negative				WMRL
Opiate Screen	Negative				WMRL

The opiate screen interpretation result includes immunoassay screen tests for both general opiates and oxycodone.

Phencyclidine Screen Negative WMRL Propoxyphene Screen Negative WMRL THC Screen Negative WMRL Tramadol Negative WMRL Creatinine WMRL 75 mg/dL Adulterant Negative WMRL See Below See Below WMRL

Decision Limits

Screen	Cutoff	Conc.	Confirm Cutoff	Conc.
Amphetamine	500	ng/mL	(GC/MS) 100	ng/mL
Barbiturates	200	ng/mL	(GC/MS) *	ng/mL
Benzodiazepines	100	ng/mL	(LC/MS/MS) *	ng/mL
Buprenorphine	5	ng/mL	(LC/MS/MS) 5	ng/mL
Cocaine	150	ng/mL	(GC/MS) 30	ng/mL
Ethanol	20	mg/dL	(GC/FID) 20	mg/dL
Fentanyl	2	ng/mL	(LC/MS/MS) 1	ng/mL
Meperidine	200	ng/mL	(LC/MS/MS) 5	ng/mL
Methadone	150	ng/mL	(GC/MS) 100	ng/mL
Opiates	300	ng/mL	(LC/MS/MS) 25	ng/mL
Oxycodone	100	ng/mL	(LC/MS/MS) 25	ng/mL
Phencyclidine	25	ng/mL	(GC/MS) 5	ng/mL
Propoxyphene	300	ng/mL	(GC/MS) 60	ng/mL
THC (Cannabis)	25	ng/mL	(GC/MS) 3	ng/mL
Tramadol	200	ng/mL	(LC/MS/MS) 5	ng/mL

LAB: L - LOW, H - HIGH, AB - ABNORMAL, C - CRITICAL,

Ordered By: KAJAL V SITWALA, MD F326000082

 $\begin{array}{c} \text{WX0000003827} \\ \text{Printed D\&T: } 9/27/2023 \ 8:42 \ AM \end{array}$

WMB-23-879 PAGE 1 OF 2



LABORATORY REPORT

Example Client, XYZ123 1234 Warde Road Ann Arbor MI 48108 **EXAMPLE, REPORT W** WX0000003827 M 07/08/1978

Test Name Result Flag Ref-Ranges Units Site

* Decision limit depends on drug detected.

Adulterant Decision Limit:
General Oxidants 200 ug/mL

The adulterant assay tests for General Oxidants, including Chromates and Nitrites. Adulterants are substances either ingested or added directly to a urine specimen to prevent the detection of drug use.

If applicable, any drug confirmation testing reported here was developed and the performance characteristics determined by Warde Medical Laboratory. This confirmation testing has not been cleared or approved by the FDA. The laboratory is regulated under CLIA as qualified to perform high-complexity testing. This test is used for patient testing purposes. It should not be regarded as investigational or for research.

Performing Site:

WMRL: Warde Medical Laboratory 300 West Textile Road Ann Arbor MI 48108 (800)876-6522

LAB: L - LOW, H - HIGH, AB - ABNORMAL, C - CRITICAL,

F326000082 Ordered By: KAJAL V SITWALA, MD WX0000003827 WX00000000002365 Printed D& T: 9/27/2023 8:42 AM

WMB-23-879 PAGE 2 OF 2