

Update Summary

Inactivate Test With Replacement	1/30/2024	MMMNG - "Myeloid Malignancies Mutation Panel by Next Gen. Sequencing" replaced by MNGS - "Myeloid Neoplasm NGS Panel"
---	-----------	---

Inactivate Test With Replacement

Effective Date	1/30/2024
-----------------------	-----------

Inactivated Test

Name	Myeloid Malignancies Mutation Panel by Next Gen. Sequencing
Code	MMMNG
Legacy Code¹	MMMNG
Interface Order Code	3618421
Notes	Inactivate ARUP MMMNG test.

Replacement Test

Name	Myeloid Neoplasm NGS Panel
Code	MNGS
CPT Code(s)	81455
Notes	Bringing NGS testing in-house

Specimen Requirements

Specimen Required	<p>Collect: Whole Blood EDTA Lavender</p> <p>Specimen Preparation: Send 1.0mL whole blood in a screw capped plastic vial.</p> <p>Minimum Volume: 0.50mL</p> <p>Transport Temperature: Refrigerated</p> <p>New York Approval Status: No</p>
--------------------------	--

Alternate Specimen	Bone Marrow
---------------------------	-------------

Rejection Criteria	Serum, plasma, tissue, buccal brush or swab, grossly hemolyzed specimens
---------------------------	--

Stability	Room Temperature: 3 days Refrigerated: 7 days Frozen: 30 days
------------------	---

Performing Information

Methodology	Targeted Next Generation Sequencing
--------------------	-------------------------------------

Reference Range	See Report
------------------------	------------

Performed Days	Monday - Friday
-----------------------	-----------------

Turnaround Time	7-14 days
------------------------	-----------

Performing Laboratory	Warde Medical Laboratory		
Interface Information			
Legacy Code¹	MNGS		
Interface Order Code	3000886		
Result Code	Name	LOINC Code	AOE/Prompt²
3000886	MNGS Report	35474-6	No