

### Update Notes

Test code **VAR – Varicella Zoster IgG Antibody**: Due to a manufacturer discontinuing their reagent, this test will be inactivated and replaced with **VARG – Varicella Zoster IgG Antibody**.

Test code **VARGM – Varicella IgG/IgM Antibodies**: Due to a manufacturer discontinuing their reagent, this test will be inactivated and replaced with **VZVGM – Varicella IgG/IgM Antibodies**.

### Update Summary

Inactivate Test With Replacement	9/11/2024	<a href="#">VAR - "Varicella Zoster IgG Antibody" replaced by VARG - "Varicella Zoster IgG Antibody"</a>
Inactivate Test With Replacement	9/11/2024	<a href="#">VARGM - "Varicella IgG/IgM Antibodies" replaced by VZVGM - "Varicella IgG/IgM Antibodies"</a>

Inactivate Test With Replacement			
<b>Effective Date</b>	9/11/2024		
Inactivated Test			
<b>Name</b>	Varicella Zoster IgG Antibody		
<b>Code</b>	VAR		
<b>Legacy Code</b>	VAR		
<b>Interface Order Code</b>	3008000		
Replacement Test			
<b>Name</b>	Varicella Zoster IgG Antibody		
<b>Code</b>	VARG		
<b>CPT Code(s)</b>	86787		
Specimen Requirements			
<b>Specimen Required</b>	<i>Collect:</i> Serum separator tube (SST) <i>Specimen Preparation:</i> Centrifuge, separate serum from cells and send 1.0 mL serum in a screw capped plastic vial. <i>Minimum Volume:</i> 0.5 mL <i>Transport Temperature:</i> Refrigerated		
<b>Alternate Specimen</b>	Plasma: Heparin, sodium and lithium; K2-EDTA		
<b>Rejection Criteria</b>	Grossly hemolyzed or lipemic samples; obvious microbial contamination; heat inactivated samples		
<b>Stability</b>	Room temperature: 3 days Refrigerated: 7 days Frozen: 1 month		
Performing Information			
<b>Methodology</b>	Chemiluminescence Immunoassay		
<b>Reference Range</b>	<1.00 S/CO Negative      >=1.00 S/CO Positive		
<b>Performed Days</b>	Sunday - Friday		
<b>Turnaround Time</b>	1 -2 days		
<b>Performing Laboratory</b>	Warde Medical Laboratory		
Interface Information			
<b>Legacy Code</b>	VARG		
<b>Interface Order Code</b>	3000904		
<b>Result Code</b>	<b>Name</b>	<b>LOINC Code</b>	<b>AOE/Prompt</b>
3000904	Varicella Zoster IgG Antibody	8047-3	No



# LABORATORY REPORT

Example Client, XYZ123  
1234 Warde Road  
Ann Arbor MI 48108

**EXAMPLE, REPORT W**  
WX0000003826 F 12/05/1988 35 Y

## Immunology

Collected: 09/06/2024 14:28 Received: 09/06/2024 14:28

<u>Test Name</u>	<u>Result</u>	<u>Flag</u>	<u>Ref-Ranges</u>	<u>Units</u>	<u>Site</u>
Varicella Zoster IgG Antibody	0.50		<1.00	S/CO	WMRL

**INTERPRETATION:**

Negative: No antibody detected.

**Reported Date:** 09/06/2024 14:28 VARG

Performing Site:

WMRL: WARDE MEDICAL LABORATORY 300 West Textile Road Ann Arbor MI 48108

LAB: L - LOW, H - HIGH, AB - ABNORMAL, C - CRITICAL, . - NOT TESTED

G506000004  
WX0000003826

Ordered By: KAJAL SITWALA, MD, PHD  
WX00000000002353

Printed D&T: 09/06/24 14:29

Kajal V. Sitwala, MD, PhD - Medical Director

Form: MM RL1

PAGE 1 OF 1

Inactivate Test With Replacement			
<b>Effective Date</b>	9/11/2024		
Inactivated Test			
<b>Name</b>	Varicella IgG/IgM Antibodies		
<b>Code</b>	VARGM		
<b>Legacy Code</b>	VARGM		
<b>Interface Order Code</b>	3017440		
Notes			
The Varicella Zoster IgM Antibody (3017400) is not changing.			
Replacement Test			
<b>Name</b>	Varicella IgG/IgM Antibodies		
<b>Code</b>	VZVGM		
<b>CPT Code(s)</b>	86787 x2		
Specimen Requirements			
<b>Specimen Required</b>	<i>Collect:</i> Serum separator tube (SST) <i>Specimen Preparation:</i> Centrifuge, separate and send 1.0 mL serum in a screw capped plastic vial. <i>Minimum Volume:</i> 0.5 mL <i>Transport Temperature:</i> Refrigerated		
<b>Alternate Specimen</b>	Plasma: Heparin, sodium and lithium; K2-EDTA		
<b>Rejection Criteria</b>	Grossly hemolyzed or lipemic samples; obvious microbial contamination; heat inactivated samples		
<b>Stability</b>	Room temperature: 3 days Refrigerated: 7 days Frozen: 1 month		
Performing Information			
<b>Methodology</b>	Chemiluminescence Immunoassay		
<b>Reference Range</b>	<1.00 S/CO Negative      >=1.00 S/CO Positive		
<b>Performed Days</b>	Sunday - Friday		
<b>Turnaround Time</b>	2 – 4 days		
<b>Performing Laboratory</b>	Warde Medical Laboratory		
Interface Information			
<b>Legacy Code</b>	VZVGM		
<b>Interface Order Code</b>	3000906		
Result Code	Name	LOINC Code	AOE/Prompt
3000904	Varicella Zoster IgG Antibody	8047-3	No
3017400 (current test)	Varicella Zoster IgM Antibody	8048-1	No