

Example Client, XYZ123 1234 Warde Road Ann Arbor MI 48108 **EXAMPLE, REPORT W** WX0000003827 M 07/08/1978 45 Y

Referral Testing					
		9/08/2023 09:46	Received:	09/08/2023	09:46
Test Name	Result	Flag Ref-Rang	ges <u>L</u>	<u>Jnits</u>	<u>Site</u>
Cortisol, Saliva, LC/MS/MS	0.35		n	ncg/dL	QCRL
 8-10 AM: 0.04-0.56 mcg/dL noon-2 PM: < OR = 0.21 mcg/dL 4-6 PM: < OR = 0.15 mcg/dL 10 PM-1 AM: < OR = 0.09 mcg/dL Please note, Salivette(R) (Order number: 51.1534) is not an acceptable saliva collection device for this assay. We will reject all white top Salivette(R) collection devices as wrong specimen type. Please contact Quest PSC for the Salivette(R) Cortisol with blue screw top (Order number: 51.1534.500) if you are still using white top Salivette(R). This test was developed and its analytical performance characteristics have been determined by Quest Diagnostics Nichols Institute San Juan Capistrano. It has not been cleared or approved by FDA. This assay has been validated pursuant to the CLIA regulations and is used for clinical purposes. Test Performed at: Quest Diagnostics Nichols Institute 33608 Ortega Highway San Juan Capistrano, CA 92675-2042 I Maramica MD, PhD, MBA 					
Performing Site: QCRL: QUEST DIAGNOSTICS REFERENCE LAB CAPISTRANO 33608 Ortega Highway San Juan Capistrano CA 92675					

Reported Date: 2023.09.08 9:46

LAB: L - LOW, H - HIGH, AB - ABNORMAL, C - CRITICAL, . - NOT TESTED