

Update Notes

Update Summary

Inactivate Test With Replacement	12/14/2021	CDHLA - "Celiac Disease Genotyping" replaced by CDHDG - "Celiac Disease HLA-DQ Genotyping"
---	------------	--

Inactivate Test With Replacement

Effective Date	12/14/2021
-----------------------	------------

Inactivated Test

Name	Celiac Disease Genotyping
Code	CDHLA
Legacy Code¹	CDHLA
Interface Order Code	3621300
Notes	

Replacement Test

Name	Celiac Disease HLA-DQ Genotyping
Code	CDHDG
CPT Code(s)	81382 x 2
Notes	

Specimen Requirements

Specimen Required	<p><i>Collect:</i> Lavender EDTA</p> <p><i>Specimen Preparation:</i> Send 3.0 mL whole blood refrigerated.</p> <p><i>Minimum Volume:</i> 1.0 mL</p> <p><i>Transport Temperature:</i> Refrigerated</p>
Alternate Specimen	Whole blood: Yellow ACD A
Rejection Criteria	Yellow ACD B. Clotted, grossly hemolyzed, or heparinized specimens.

Stability	Room temperature: 72 hours Refrigerated: 7 days Frozen: Unacceptable		
Performing Information			
Methodology	Polymerase Chain Reaction/Massively Parallel Sequencing Polymerase Chain Reaction/Sequence-Specific Oligonucleotide Probe		
Reference Range	See report		
Performed Days	Monday - Friday		
Turnaround Time	10 - 18 days		
Performing Laboratory	ARUP Reference Laboratory		
Interface Information			
Legacy Code¹	CDHDG		
Interface Order Code	3600221		
Result Code	Name	LOINC Code	AOE/Prompt²
3600222	HLA-DQA1, Allele 1	59019-0	No
3600223	HLA-DQA1, Allele 2	59019-0	No
3600224	HLA-DQB1, Allele 1	57299-0	No
3600225	HLA-DQB1, Allele 2	57299-0	No
3600226	Celiac HLA Interpretation	48767-8	No