



LABORATORY REPORT

QC ACCOUNT (WARDE)
300 W. TEXTILE
ANN ARBOR MI 48108

EXAMPLE, REPORT

WX0000000237 F 12/05/1988 36 Y

Referral Testing

Collected: 09/18/2025 09:34

Received: 09/18/2025 09:34

Test Name	Result	Flag	Ref-Ranges	Units	Site
Soluble Transferrin Receptor	2.5	H	0.90-2.01	mg/L	ARRL

INTERPRETIVE INFORMATION: Soluble Transferrin Receptor

The Beckman Coulter Access sTfR immunoassay is intended as an aid in the diagnosis of iron deficiency anemia, especially in patients with chronic disease. In adult patients with anemia, an sTfR result greater than or equal to 1.55 mg/L is 86 percent sensitive and 49 percent specific for the presence of iron deficiency anemia, alone or in combination with anemia of chronic disease. The sTfR assay is not intended to be used in isolation; results should be interpreted in conjunction with the patient's clinical presentation and other diagnostic tests, such as other indicators of iron status (refer to table below).

Analyte	Tests for Changes in:	Iron Def. Anemia	Anemia of Chronic Disease	Combined Iron Def. and anemia of Chronic Dz
Ferritin	Iron Stores	Low	High	Normal or High
TIBC	Iron Status	High	Low	Normal or High
Serum Fe	Iron Status	Low	Low	Low
sTfR	Iron Status	High	Normal	High

Performed By: ARUP Laboratories
500 Chipeta Way

Salt Lake City, UT 84108

Laboratory Director: Jonathan R. Genzen, MD, PhD

CLIA Number: 46D0523979

Reported Date: 09/18/2025 09:34 STR

Performing Site:

ARRL: ARUP REFERENCE LAB 500 Chipeta Way Salt Lake City UT 841081221

LAB: L - LOW, H - HIGH, AB - ABNORMAL, C - CRITICAL, . - NOT TESTED

E118000007
WX0000000237

Printed D&T: 09/18/25 09:34

Ordered By: CLIENT CLIENT
WX00000000000336

Kajal V. Sitwala, MD, PhD - Medical Director

Form: MM RL1

PAGE 1 OF 1